

microlife®

A partner for people. For life.



▶ Product video



Atrial Fibrillation detection with

AFIBsens

The first step in stroke prevention.

BP B3 AFIB Blood Pressure Monitor

- Atrial fibrillation detection for stroke prevention
- MAM: 3 measurements in 1 click & average result
- Comfortable cuff for most arm sizes (22–42 cm)
- Cuff fit check
- Date and time

Key features



High accuracy in Atrial Fibrillation (AF) detection for stroke prevention



IHB for irregular heartbeat detection



Gentle+ ensures optimal speed and pressure control for comfortable measurements



Easy to read classification of blood pressure result



Extra large display



99 data memory and average-display of all stored readings



Cuff fit check



MAM: Microlife Average Mode

- 3 measurements in 1 click and average result for superior reliability
- Reliable hypertension diagnosis and therapy control



Soft Cuff: 22-42 cm

- Conical shape for better fit and optimal accuracy
- Washable cuff cover for longlasting, hygienic use



Clinically Validated

Highest measuring accuracy with clinical A/A grade according to the British and Irish Hypertension Society (BIHS).

12+ years

Validated and recommended for children. Also validated and recommended for patients with diabetes, renal disease, pregnancy / pre-eclampsia and obesity.



5 year guarantee