

# microlife®

A partner for people. For life.

Atrial Fibrillation detection with

# AFIBsens

The first step in stroke prevention.



reddot winner 2021

Bluetooth®

## BP A7 Touch BT Blood Pressure Monitor

- Atrial fibrillation detection for stroke prevention
- Touch screen technology
- Easy data analysis with PC software included and app for smartphone



Product video



### Key features



High accuracy in Atrial Fibrillation (AF) detection for stroke prevention



Microlife Average Mode (MAM): 3 measurements in 1 click



Gentle+ ensures optimal speed and pressure control for comfortable measurements



Easy to read classification of blood pressure result



Large touch screen



Comfortable, preformed cuff for most arm sizes (22-42 cm)



Two users with guest mode and 2 x 99 data memory



#### Microlife Connected Health+ App

- Bluetooth® 4.0 for automatic data transfer after each measurement
- Compatibility: Android™ OS 5 or newer / iPhone iOS 10.0 or newer
- More on [www.microlife.com/connect](http://www.microlife.com/connect)



#### MAM: Microlife Average Mode

- 3 measurements in 1 click and average result for superior reliability
- Reliable hypertension diagnosis and therapy control



#### Clinically Validated

Highest measuring accuracy with clinical A/A grade according to the British and Irish Hypertension Society (BIHS).



Validated and recommended for children. Also validated and recommended for patients with diabetes, renal disease, pregnancy / pre-eclampsia and obesity.



The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Microlife Corp. is under license. Other trademarks and trade names are those of their respective owners.



5 year  
guarantee