

microlife®

ASTHMA MONITOR

PEF & FEV₁
measurement



PF100

Digital peak flow meter

- Suitable for children & adults
- High accuracy
- Asthma Analyser Software for the management of the patient's respiratory diary

Key features



High accuracy
Peak Expiratory Flow (PEF) & FEV₁
measurement



240 data memory
Stored with date and time



Traffic light system
For an immediate assessment
of peak expiratory flow



PC connection with Asthma Analyser
Software (free software download),
transfer cable included



Rotated vane wheel method



Suitable for children and adults



Professional carrying bag

